



ECKERMANN

Triplo®

DENTAL IMPLANT SYSTEM



Implantes

Triplo®

Tapón de cierre incluido

- Ø 4,0 mm. REF.
- ‡ 8,0 mm. IKT 354008
- ‡ 10,0 mm. IKT 354010
- ‡ 11,5 mm. IKT 354011
- ‡ 13,0 mm. IKT 354013
- ‡ 15,0 mm. IKT 354015
- ‡ 17,0 mm. IKT 354017

- Ø 4,5 mm.
- ‡ 6,0 mm. IKT 354506
- ‡ 8,0 mm. IKT 354508
- ‡ 10,0 mm. IKT 354510
- ‡ 11,5 mm. IKT 354511
- ‡ 13,0 mm. IKT 354513
- ‡ 15,0 mm. IKT 354515
- ‡ 17,0 mm. IKT 354517

- Ø 5,0 mm.
- ‡ 6,0 mm. IKT 355006
- ‡ 8,0 mm. IKT 355008
- ‡ 10,0 mm. IKT 355010
- ‡ 11,5 mm. IKT 355011
- ‡ 13,0 mm. IKT 355013
- ‡ 15,0 mm. IKT 355015
- ‡ 17,0 mm. IKT 355017

Triplo® perio

Tapón de cierre incluido

- Ø 4,0 mm. REF.
- ‡ 8,0 mm. IPT 354008
- ‡ 10,0 mm. IPT 354010
- ‡ 11,5 mm. IPT 354011
- ‡ 13,0 mm. IPT 354013
- ‡ 15,0 mm. IPT 354015
- ‡ 17,0 mm. IPT 354017

- Ø 4,5 mm.
- ‡ 6,0 mm. IPT 354506
- ‡ 8,0 mm. IPT 354508
- ‡ 10,0 mm. IPT 354510
- ‡ 11,5 mm. IPT 354511
- ‡ 13,0 mm. IPT 354513
- ‡ 15,0 mm. IPT 354515
- ‡ 17,0 mm. IPT 354517

- Ø 5,0 mm.
- ‡ 6,0 mm. IPT 355006
- ‡ 8,0 mm. IPT 355008
- ‡ 10,0 mm. IPT 355010
- ‡ 11,5 mm. IPT 355011
- ‡ 13,0 mm. IPT 355013
- ‡ 15,0 mm. IPT 355015
- ‡ 17,0 mm. IPT 355017

Instrumental

Caja Master Box

- REF. CKS

Caja Starter Box

- REF. CKSX

Caja Basic Box

- REF. CKSB

Caja Inox Box

- REF. CKO

Bisturí circular

- REF. MUF 500
- Ø 5,0 mm. MUF 550

Fresa lanceolada portabisturí

- REF. FLM 200
- 2,0 mm.

Fresa lanceolada

- REF. EX 408
- REF. EX 409
- Ø corta
- larga

Fresa de avellanado

- REF. EX 407

Fresa piloto

- REF. EX 600

Fresa helicoidal de acabado

- REF. FSH 25
- Ø 2,5 mm.
- REF. FSH 30
- Ø 3,0 mm.
- REF. FSH 35
- Ø 3,5 mm.
- REF. FSH 40
- Ø 4,0 mm.
- REF. FSH 45
- Ø 4,5 mm.
- REF. FSH 50
- Ø 5,0 mm.

Fresa helicoidal de acabado corta

- REF. FSH 25 C
- Ø 2,5 mm.
- REF. FSH 30 C
- Ø 3,0 mm.
- REF. FSH 35 C
- Ø 3,5 mm.
- REF. FSH 40 C
- Ø 4,0 mm.
- REF. FSH 45 C
- Ø 4,5 mm.
- REF. FSH 50 C
- Ø 5,0 mm.

Topes de fresado

- REF. TFSH 6
- ‡ 6,0 mm.
- REF. TFSH 8
- ‡ 8,0 mm.
- REF. TFSH 10
- ‡ 10,0 mm.
- REF. TFSH 11
- ‡ 11,5 mm.
- REF. TFSH 13
- ‡ 13,0 mm.
- REF. TFSH 15
- ‡ 15,0 mm.
- REF. TFSH 17
- ‡ 17,0 mm.

Fresa helicoidal de acabado Golden

- REF. FIS 25
- Ø 2,5 mm.
- REF. FIS 30
- Ø 3,0 mm.
- REF. FIS 35
- Ø 3,5 mm.
- REF. FIS 40
- Ø 4,0 mm.
- REF. FIS 45
- Ø 4,5 mm.
- REF. FIS 50
- Ø 5,0 mm.

Topes de fresado Golden

- REF. TFSH 6
- ‡ 6,0 mm.
- REF. TFSH 8
- ‡ 8,0 mm.
- REF. TFSH 10
- ‡ 10,0 mm.
- REF. TFSH 11
- ‡ 11,5 mm.
- REF. TFSH 13
- ‡ 13,0 mm.
- REF. TFSH 15
- ‡ 15,0 mm.
- REF. TFSH 17
- ‡ 17,0 mm.

Prolongador de fresa

- REF. EX 650

Paralelizador

- REF. PL 15

Medidor de profundidad

- REF. MP 25

Iniciador de rosca

- REF. MRS 2010
- Ø 2,0 mm.
- REF. MRS 2510
- Ø 2,5 mm.
- REF. MRS 3010
- Ø 3,0 mm.
- REF. MRS 3510
- Ø 3,5 mm.
- REF. MRS 4010
- Ø 4,0 mm.
- REF. MRS 4510
- Ø 4,5 mm.
- REF. MRS 5010
- Ø 5,0 mm.

Llave rosca implantes

- REF. LI IKT
- REF. LI IKT L
- Ø corta
- larga

Llave rosca implantes con contra-ángulo

- REF. LC IKT

Carraca

- REF. CR

Carraca dinamométrica

- REF. CRT 35
- REF. CRT 45
- 35 Nw
- 45 Nw

Destornillador hexagonal

- REF. DH 17 T
- REF. DH 23 T
- Ø corto
- largo

Llave portamechas universal

- REF. LMU

Mecha hexagonal

- REF. MH

Mecha destornillador cruceta

- REF. MC 18

Instrumento portamechas

- REF. DM

Conformación de tejidos

Tornillo de cierre

- REF. TC IKT

Moldeador gingival

- REF. MG 5515 IKT
- REF. MG 5530 IKT
- REF. MG 5550 IKT
- REF. MG 5570 IKT
- ‡ 1,5 mm.
- ‡ 3,0 mm.
- ‡ 5,0 mm.
- ‡ 7,0 mm.
- Ø 5,5 mm.

Moldeador gingival periodontal

- REF. MG 5515 IPT*
- REF. MG 5530 IPT
- REF. MG 5550 IPT
- REF. MG 5570 IPT
- ‡ 1,5 mm.
- ‡ 3,0 mm.
- ‡ 5,0 mm.
- ‡ 7,0 mm.
- Ø 5,5 mm.

Toma de impresión

Transfer de impresión directa + tornillo

- REF. TR IKT P
- Ø 5,5 mm.

Transfer de impresión indirecta + tornillo

- REF. TR IKT A
- Ø 5,5 mm.

Réplica implante Triplo / Triplo perio

- REF. RPT 40
- REF. RPT 45
- REF. RPT 50
- Ø 4,0 mm.
- Ø 4,5 mm.
- Ø 5,0 mm.

Prótesis

Tornillo oclusal

- REF. THX 1813

Tornillo de asiento

- REF. TCM 40 IKT
- REF. TCM 60 IKT
- ‡ 4,0 mm.
- ‡ 6,0 mm.

Pilar mecanizado ConeLock

- REF. PM 4516 IKT
- REF. PM 6516 IKT
- REF. PM 9516 IKT
- Ø 5,0 mm.
- ‡ Hombro 1,6 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.
- ‡ 9,0 mm.

Pilar mecanizado periodontal

- REF. PM 4525 IKT
- REF. PM 6525 IKT
- REF. PM 9525 IKT
- Ø 5,0 mm.
- ‡ Hombro 2,5 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.
- ‡ 9,0 mm.

Pilar mecanizado periodontal

- REF. PM 4516 IPT
- REF. PM 6516 IPT
- REF. PM 9516 IPT
- Ø 5,0 mm.
- ‡ Hombro 1,6 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.
- ‡ 9,0 mm.

Pilar mecanizado macizo

- REF. PMC 5040 IKT
- REF. PMC 5060 IKT
- REF. PMM 5040 IKT
- REF. PMM 5060 IKT
- Ø 5,0 mm.
- ‡ Hombro 1,6 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.
- ‡ Hombro 2,5 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.

Pilar mecanizado angulado

- REF. PMA 20 IKT*
- Ø 5,0 mm.
- ‡ ANGLADÓN 20°
- ‡ 9,0 mm.

Pilar calcinable

- REF. CA IKT 50 A
- REF. CA IKT 50 C*
- Ø 5,0 mm.
- ‡ 9,0 mm.
- Ø 5,0 mm.
- ‡ 9,0 mm.

Asiento calcinable

- REF. AC

Cofia calcinable

- REF. CA PM 40
- REF. CA PM
- Ø 5,5 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.

Asiento metálico

- REF. AM

Cofia metálica

- REF. CM PM 40
- REF. CM PM
- Ø 5,5 mm.
- ‡ 4,0 mm.
- ‡ 6,0 mm.

Pilar Locator

- REF. PL IKT 50*
- Ø 4,5 mm.
- ‡ 5,0 mm.

Convertor a hexágono

- REF. C15 IKT IH
- REF. C27 IKT IH
- REF. C42 IKT IH
- ‡ 1,5 mm.
- ‡ 2,7 mm.
- ‡ 4,2 mm.



DENTAL IMPLANT SYSTEM

Avda. Extremadura, 4 03300 Orihuela - Alicante (España)
Tlf. 902 306 464 Fax. 966 741 229 E-Mail: ecklab@eckermann.es
www.eckermann.es

